

Supplementary material has been published as submitted. It has not been copyedited, or typeset by Acta Dermato-Venereologica

Table SI.

CD4 ⁺ T cell isolation kit	Miltenyi Biotec, Bergisch Gladbach, North Rhine-Westphalia, Germany
Leukoreduction chambers	NHS Blood and Transplant, Cambridge, UK
PBMCs	Cellero, Lowell, MA, USA
anti-OX40L (amlitelimab)	Kymab, Cambridge, Cambridgeshire, UK
anti-OX40 hIgG1 (A)	In-house generated
anti-OX40 hIgG1 (B)	In-house generated
hIgG4(PE) isotype control	In-house generated
hIgG1 isotype control	In-house generated
anti-hCD4	Clone: SK3; Biolegend, San Diego, CA, USA
anti-hCD25	Clone: M-A251; Biolegend, San Diego, CA, USA
anti-hFOXP3	Clone: 206D; Biolegend, San Diego, CA, USA
anti-hCD127	Clone: A019D5; Biolegend, San Diego, CA, USA
NK cell isolation kit	Miltenyi Biotec, Bergisch Gladbach, North Rhine-Westphalia, Germany
anti-OX40L hIgG1 (C)	Abnova, Taipei, Taiwan
anti-OX40L hIgG1 (D)	In-house generated
DELFIATM cell cytotoxicity kit	PerkinElmer, Shelton, CT, USA

Ig, immunoglobulin; HEK, human embryonic kidney, NK, natural killer, PBMC, peripheral blood mononuclear cell.